

802DY Series

aluminum high strength, lightweight flip-back-arm wheelchairs

- Medicare: K0004
- Folding high strength aluminum frame in black powder coat finish
- Padded flip-back desk armrests with plastic side panels
- Swing-away detachable footrests with heel strap
- 24" x 1" mag wheels with polyurethane tires
- 8" x 1" caster wheels with polyurethane tires
- Adjustable length anti-slip footplates.
- High-grade flame retardant, breathable nylon upholstery

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	weight w/o riggings (lb.)	weight Capacity (lb.)
802DY-18	18	16	20	17	37	25	12	30	250
802DY-16	16	16	20	16	37	23	12	30	250

Model No.	Description
802DY-18	18", aluminum, high strength lightweight, flip-back-arms, swing-away footrests
802DY-18-E	18", aluminum, high strength lightweight, flip-back-arms, elevating legrests
802DY-16	16", aluminum, high strength lightweight, flip-back-arms, swing-away footrests
802DY-16-E	16", aluminum, high strength lightweight, flip-back-arms, elevating legrests

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Elevating Legrests Anti-Tippers Brake Extension IV Pod Attachment Oxygen Tank Holder Seat Belt Leg Straps
- Quick Release Axle Spoke Wheels Run Flat Inserts Desk Tray Carrying Pouches •

