



KM-3520F 2-in-1 Model

Lightweight aluminum high-strength w/ dinner tray/side panel

- Folding 6061 T-6 aircraft-grade aluminum frame with powder coat finish in black colors
- 2-in-1 dinner tray/side panel
- Adjustable table depth
- Folding Backrests
- Padded flip-back armrests
- Swing-away footrests with velcro mount leg strap
- 24" x 1" mag wheels with polyurethane tires
- 7" x 1" caster wheels with rubber tires
- Adjustable length anti-slip footplates.
- High-grade flame retardant, breathable nylon upholstery

2-in-1 Tray-Side Panels



Folding Backrest



Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM-3520-18	18	16	20	16	38	28	29	250
KM-3520-16	16	16	20	16	38	26	29	250

Model No.	Description
KM-3520-18	18", aluminum, high strength lightweight, flip-back-arms, swing-away footrests
KM-3520-18-E	18", aluminum, high strength lightweight, flip-back-arms, elevating legrests
KM-3520-16	16", aluminum, high strength lightweight, flip-back-arms, swing-away footrests
KM-3520-16-E	16", aluminum, high strength lightweight, flip-back-arms, elevating legrests

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Elevating Legrests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Heel Loops • Leg Straps • Spoke Wheels • Pneumatic Tires • Run Flat Inserts • Carrying Pouches •

 **KARMA HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335