

KM-5000 Series

Reclining, lightweight aluminum detachable-arm wheelchair/transport



- Folding 6061 T-6 Aircraft-grade aluminum frame in black powder coat finish
- Medicare Code: K0004+E1226+E0966
- Backrest Reclines from 90°~163°
- Full length/padded detachable armrests with side panels
- Swing-away/detachable elevating footrests
- Adjustable length leg supports & footplates
- Adjustable-height head pillow
- 24" x 1" Mag wheels with flat free polyurethane tires
- Flat free front polyurethane casters
- Adjustable toggle wheel locks
- Includes anti-tippers
- High-grade flame retardant, breathable nylon upholstery
- Limited lifetime warranty on frames and cross-braces, one year on parts
- One of only two lightweight recliners available in the USA
- Also Available in Small Transporter Rear Wheel option, (KM-5000TP), please call for more info

(*Dimensions shown are for Large 24" Wheel Model)

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM-5000-16	16	18	20	33 (19+15)	51	25	36	250
KM-5000-18	18	18	20	33 (19+15)	51	27	36	250
KM-5000-20	20	18	20	33 (19+15)	51	29	38	300
KM-5000-22	22	18	20	33 (19+15)	51	31	38	350

** Price & specifications are subject to change without notice



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335