

KM-802 Series

*Aluminum high strength, lightweight
flip-back armrest*



- K0004 & K0005 Medicare Code
- 6061 T-6 Aircraft-grade aluminum frame
- Padded full sized/flip-back armrests with side panels
- Swing-away, detachable, height adjustable footrests with leg strap
- Dual axle for adjusting seat height to 18 or 19 in. (fixed axle model only)
- 24" x 1" spoke wheels with high profile polyurethane tires (flat free)
- Optional: Quick-release rear wheels
- 8" x 1" polyurethane flat free front casters
- High-grade flame retardant, breathable nylon upholstery
- Easy tool-less removable footplates
- Black powder coat finish
- Easy push manual brakes
- Weight capacity: 250 lbs
- Lifetime Limited Warranty on Frame

SPECIFICATIONS (subject to change without notice) :

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM802F16	16	17	9	18-19	12/18	38	23	12	42	30	250
KM802F18	18	17	9	18-19	12/18	38	25	12	42	30	250
KM802Q16	16	17	9	19	12/18	38	23	12	42	30	250
KM802Q18	18	17	9	19	12/18	38	25	12	42	30	250

OPTIONIONAL ACCESSORIES:

- Anti-Tippers - Brake Extension - IV Pod Attachment - Oxygen Tank Holder - Seat Belt - Leg Straps
 - Headrest - Backrest Extension - Carrying Pouches - Frog Leg Suspensions - UV Reflecting Canopy- Elevating Legrest



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335