











KM-BT10

Adaptive Bariatric Manual Wheelchair

- Foldable and robust frame design, easy to stow and carry
- Smooth hand-rims for easy self-propelling
- Fully detachable & flip-back armrests for ultimate convenience, also height adjustable
- Ergonomic arm-pads which provide a proper position for arm placement
- Press-down easy-tab for fast & safe detachment of swing-away footrests (enables patient to easily remove footrests without assistance)
- Attendant Brake & Manual Toggle Brake
- Rear wheel Center of Gravity Adjustment and Height adjustable.
- Exclusive lateral support bar which provides enhanced stability and energy-saving when propelled by attendant
- Angle/Height Adjustable Backrest w/ Tension Adjustable Strap
- Angle/Height/Depth Adjustable Footplate
- 24" x 2 1/8" Flat Free Rear Wheels
- 7" x 1 3/4" Flat Free Front Casters
- Frame Color: Diamond Black
- Limited Lifetime Warranty on Frame
- Chair Only: 63-70 lbs (w/o Footrests)
- Weight Capacity For All Sizes: 550 lbs



items	F/R wheels	Seat width	Seat depth	Overall	Folded	Arm height	Overall	Seat to Floor	Back height	Overall	Weight	Max. User Weight
Width												
BT10-22W	7 / 24	22	18 / 20	31	19	10½ - 13½	47	15.5-19.5	16-20	34 - 38	63 lbs	550 lbs.
BT10-24W	7 / 24	24	18 / 20	33	19	10½ - 13½	47	15.5-19.5	16-20	34 - 38	65 lbs	550 lbs.
BT10-26W	7 / 24	26	18 / 20	35	19	10 ½ - 13½	47	15.5-19.5	16-20	34 - 38	67 lbs	550 lbs.
BT10-26W	7 / 24	26	22 / 24	35	19	10 ½ - 13½	47	15.5-19.5	16-20	34 - 38	67 lbs	550 lbs.
BT10-28W	7 / 24	28	22 / 24	37	20	10 ½ - 13½	47	15.5-19.5	16-20	34 - 38	68 lbs	550 lbs.
BT10-30W	7 / 24	30	22 / 24	39	20	10 ½ - 13½	47	15.5-19.5	16-20	34 - 38	70 lbs	550 lbs.

** Dimensions are in inches **

Companion push bar is designed to assist any caretaker for greater maneuverability.

Flip back armrests allows for easier and safe transferring in and out of the BT-10. Height adjustable armrests are available.

Multiple adjustments such as Rear wheel positioning, seat to floor height, and back angle make BT-10 a perfect fit for anyone.



Our unique frame design is multi-adjustable and has a maximum weight capacity of 550 lbs.



Fully customizable seat allows for greater seat width, depth, and overall better positioning for any patient



Multi adjustable footplate ensures comfort and safety



Unique frame design allows for easy and compact folding



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry, CA 91748-1418
Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335