



KM-2020 Series

Aluminum high strength, lightweight flip-back-arm wheelchair

- Folding 6061 T-6 Aircraft-grade aluminum frame with powder coat finish in optional colors
- Full padded flip-back armrests with side panels
- Swing-away detachable footrests with Velcro mount leg strap
- 14"x 1" Rear Pneumatic Poly-urethane tires
- Ergonomic push grips
- Extra long adjustable toggle wheel locks
- Adjustable length anti-slip footplates
- High-grade flame retardant, breathable nylon upholstery w/ rear pouch
- 8"x 1" Front casters

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM-2020-18	18	16	20	17	36	24 ½	24	250

Frame Color Options: **Black**

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

• Elevating Leg rests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Heel Loops • Leg Straps • Quick Release Axle • Spoke Wheels • Pneumatic Tires •

