

KM-8520F-20W-HA

Aluminum high strength, lightweight flip-back-arm wheelchairs



- Folding 6061 T-6 Aircraft-grade aluminum frame with powder coat finish
- Padded & contoured height adjustable flip-back armrests with side panels
- With Folding Backrests
- Swing-away detachable footrests with Velcro mount leg strap
- 8" x 2" Heavy-duty flat-free polyurethane casters
- 24" x 1 3/4" Heavy-duty flat-free spoke wheels with polyurethane tires
- Adjustable toggle wheel locks.
- Adjustable length anti-slip footplates.
- High-grade flame retardant, breathable nylon upholstery

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM-8520W20HA	20	18	9-11	20	16	36	28	14	41	35	350

Model No. Description

- KM-8520W20HA 20", aluminum, high strength lightweight, flip-back-arms, swing-away footrests
- KM-8520W20HA-E 20", aluminum, high strength lightweight, flip-back-arms, elevating footrests

Frame Color Options: **Black**

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Elevating Footrests • Anti-Tippers • Brake Extension • Oxygen Tank Holder • Seat Belt • Leg Straps •



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335