



# KN-700T

*Standard steel detachable armrest wheelchair*

- Medicare code K0001
- Detachable/Padded desk length armrests with alloy side guard
- Swing-away/Height adjustable footrests with heel loop
- Dual-Axle/Seat height adjustable frame
- Push to lock wheel brakes
- 24" Mag wheels with flat free tires
- Extra-durable forks
- 8" Flat free front casters
- High-grade/Padded flame retardant vinyl upholstery
- Carry pocket on backrest
- Seat width: 16", 18"
- Frame Color: Silver vein
- Weight Capacity: 250 lbs

**SPECIFICATIONS (subject to change without notice) :**

	seat width	seat depth	armrest height	seat height	back height	overall height	overall open width	folded width	overall length	weight w/o riggings	weight capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
KN-700NT	16	16	9	17 or 19	16	36	24	11	31	39	250
KN-700T	18	16	9	17 or 19	16	36	26	11	31	39	250

**Optional Accessories:**

- Elevating Legrests - Anti-Tippers - Brake Extension - IV Pod Attachment - Oxygen Tank Holder - Seat Belt - Leg Straps
- Carry Pouch - Headrest - Backrest Extension



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335