



KN-800T

Standard fixed-arm wheelchairs

- Medicare Code K0001
- Fixed/Padded full sized armrests with alloy side guard
- Swing-away/Height adjustable footrests with heel loop
- 24" Mag wheels with flat free tires
- Steel extra-durable front caster forks
- 8" Flat free front caster
- Push to lock wheel brakes
- High-grade/Padded flame retardant vinyl upholstery
- Carry pocket on backrest
- Frame color: Silver vein
- Seat width: 16", 18"
- Weight capacity: 250 lbs

SPECIFICATIONS (subject to change without notice) :

	seat width	seat depth	armrest height	seat height	back height	overall height	overall open width	folded width	overall length	weight w/o riggings	weight capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
KN-800NT	16	16	11	19	16	36	23	11	31	37	250
KN-800T	18	16	11	19	16	36	25	11	31	37	250

Optional Accessories:

- Elevating Footrests - Anti-Tippers - Brake Extension - IV Pod Attachment - Oxygen Tank Holder - Seat Belt - Leg Straps -
- Carry Pouch - Headrest - Backrest Extension



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry. CA 91748-1418
 Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335