



KN-880 Series

reclining, steel, detachable arm wheelchair

- Medicare Code: K0003+K0195+E1226+E096
- Folding chrome plated steel frame
- 180° full reclining back with anti-tippers
- Padded detachable full armrests
- Detachable elevating footrests
- 24" x 1" flat-free mag wheels
- 8" x 1" flat-free casters
- Adjustable length leg supports & footplate
- High-grade flame retardant vinyl upholstery

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity riggings (lb.)
KN-880NE	16	16	20	29 (+)	49	23	50	250
KN-880E	18	16	20	29 (+)	49	25	50	250
KN-880WE	20	16	20	29 (+)	49	27	52	300

Model No.	Description
KN-880NE	16", 180° reclining back, detachable full arms, elevating legrests
KN-880E	18", 180° reclining back, detachable full arms, elevating legrests
KN-880WE	20", 180° reclining back, detachable full arms, elevating legrests

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Seat Belt • Leg Straps • Pneumatic Tires • Run Flat Inserts • Desk Tray • Carrying Pouches •



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335