



KN-920W Series

Deluxe Steel Manual Wheelchair, 20" Wide Seat

- *HCPCS - K0007*
- *High quality steel construction*
- *Detachable padded desk armrests*
- *Swing-away footrests composite footplates w/ heel loop*
- *24" x1" Mag wheels w/ solid tires*
- *7"x1" Front casters*
- *Adjustable Seat Height*
- *Double padded vinyl upholstery*
- *300 lbs weight capacity, single cross brace*
- *Frame: Chrome*

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	Weight w/o riggings (lb.)	weight capacity (lb.)
KN-920W	20"	16"	18" - 19"	18"	36"	28"	41 lbs	300
KN-920WE	20"	16"	18" - 19"	18"	36"	28"	41 lbs	300

** Price & specifications are subject to change without notice.

- Elevating Largest • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Leg Straps
- Pneumatic Tires



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335