

KN-922W

*Deluxe Steel Manual Wheelchair with
22" Wide Seat*



- *HCPCS - K0007*
- *High quality steel construction*
- *Detachable desk armrests*
- *Swing-away footrests w/ heel loop*
- *24" x 1" Mag wheels w/ solid tires & steel hand-rim*
- *8"x1.5" Solid Front Casters*
- *Adjustable seat height*
- *Double padded vinyl upholstery*
- *Frame: Black w/ silver vein*
- *Weighs 49 lbs*
- *350 lbs Weight capacity, dual cross brace*

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	Weight riggings (lb.)	weight capacity (lb.)
KN-922W	22"	18"	18" - 19"	18"	36"	30"	49	350
KN-922WE	22"	18"	18" - 19"	18"	36"	30"	49	350

** Price & specifications are subject to change without notice.

- Elevating Legrest • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Leg Straps
- Pneumatic Tires



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335