

KN-924W

*Deluxe Steel Manual Wheelchair with
24" Wide Seat*



- *HCPCS – K0007*
- *High quality steel construction*
- *Detachable desk armrests*
- *Swing-away footrests w/ heel loop*
- *24" x 1" Mag wheels w/ solid tires & steel hand-rim*
- *8"x1.5" Solid Front Casters*
- *Adjustable seat height*
- *Double padded vinyl upholstery*
- *Frame: Chrome*
- *Weighs 57 lbs*
- *400 lbs Weight capacity, dual cross brace*

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	Weight w/o riggings (lb.)	weight capacity (lb.)
KN-924W	24"	18"	20"	16"	36"	32"	56	400
KN-924WE	24"	18"	20"	16"	36"	32"	56	400

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335