



LT-2019

*aluminum lightweight
fixed arm companion chairs*

- Folding aluminum frame with powder coat finish in optional colors
- Padded fixed armrests
- Swingaway detachable footrests with heel loops
- Foldable backrests for easy storage
- 8" x 1" caster wheels with rubber tires
- Rear caster toggle wheel locks
- Adjustable length footplates
- High-grade flame retardant, breathable nylon upholstery
- Frame colors: Black, Blue or, Burgundy

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
LT-2017	17	16	19 1/2	19 (11+8)	40	21 ¼	19	250
LT-2019	19	16	19 1/2	19 (11+8)	40	23 ¼	19	250

Model No.	Description
LT-20-17	17", fixed arms, swing-away footrests
LT-20-19	19", fixed arms, swing-away footrests

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

- Elevating Legrests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Heel Loops • Leg Straps • Quick Release Axle • Spoke Wheels • Pneumatic Tires • Run Flat Inserts • Desk Tray • Carrying Pouches •



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335