

## MVP-502TP Ergonomic V-Seat Recliner

Patented Ergonomic/Anti-Sliding "V-Seating" prevents sliding, sores and relieves overall pressure on your body! backrest reclines smoothly from 93-160 degrees.

Rigidifying/Folding push-bar provides and strengthens

back support  $\rightarrow$ 



Ergonomic Armrest allows the user to adjust the armrest p

position when reclined  $\rightarrow$ 




- "Traveling-Axle" Leg-rest automatically stretches smoothly to reduce pressure/compression on abdomen.
- 3D Elevating Leg-rest support Ergonomically Designed to allow greater micro-adjustment: up/down,
- left/right or inward/outward  $\rightarrow$



- Extending handbrakes equipped!
- Strap-Tension adjustable backrest & seat for ultimate comfort adjustments and stabilization for each user!
- Ergonomic Upholstery; Also pressure relieving 1.25" to 3" thick breathable material, offers superb comfort throughout the body
- Folds Ultra Compact & Weighs only ~ 44lbs w/o all fixtures  $\rightarrow$









## KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418 Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335