



R-4100 Series

Aluminum 4-wheeled Rollators

- Folding/Lightweight aluminum Rollator w/powder coat finish in burgundy or blue
- Adjustable height, Hand-grips with locking loop brakes on (4) 6" casters
- Padded seat and curved/folding backrest for added comfort + easy to fold and store
- Easy-remove under seat storage basket with handle
- Easier to maneuver outdoors and indoors than standard Rollators with 6" casters
- Weight w/ Basket is 11 lbs
- Weight Capacity is 250 lbs

| | Seat Width (in.) | seat pad depth (in.) | Seat pad width (in.) | seat height (in.) | backrest height (in.) | overall height (in.) | overall open width (in.) | weight (lb.) | weight capacity (lb.) |
|--------|------------------|----------------------|----------------------|-------------------|-----------------------|----------------------|--------------------------|--------------|-----------------------|
| R-4100 | 18 | 12 | 12 | 18.5 | 9.5 | 28~34 | 21-23.5 | 11 | 250 |

** Price & specifications are subject to change without notice.



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335