



## R-4600 Series

### *Aluminum 4-wheeled Rollators*

- Folding/Lightweight aluminum Rollator w/powder coat finish in burgundy or blue
- Adjustable height, Hand-grips with locking loop brakes on (4) 6" casters
- Padded seat and curved/folding backrest for added comfort + easy to fold and store
- Easy-remove under seat storage basket with handle
- Easier to maneuver outdoors and indoors than standard Rollators with 6" casters
- Weight w/ Basket is 11 lbs
- Weight Capacity is 250 lbs

	Seat Width (in.)	seat pad depth (in.)	Seat pad width (in.)	seat height (in.)	backrest height (in.)	overall height (in.)	overall open width (in.)	weight (lb.)	weight capacity (lb.)
R-4600	18	14	13	22	12	24~38	24	12	250

\*\* Price & specifications are subject to change without notice.



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335

Filename: R-4600  
Directory: C:\Users\David\Documents  
Template: C:\Users\David\AppData\Roaming\Microsoft\Templates\Normal.dotm  
Title: E-800 Series  
Subject:  
Author: Uriel  
Keywords:  
Comments:  
Creation Date: 2/2/2011 2:40:00 PM  
Change Number: 5  
Last Saved On: 8/2/2012 2:53:00 PM  
Last Saved By: David  
Total Editing Time: 16 Minutes  
Last Printed On: 8/2/2012 2:53:00 PM  
As of Last Complete Printing  
Number of Pages: 1  
Number of Words: 141 (approx.)  
Number of Characters: 808 (approx.)