

# R-4608 Series

*Aluminum 4-wheeled Rollators w/8" wheels*



- Folding/Lightweight aluminum Rollator w/powder coat finish in burgundy or blue
- Adjustable height, Hand-grips with locking loop brakes on (4)8" casters
- Padded large seat and curved/folding backrest for added comfort + easy to fold and store
- Easy-remove under seat storage basket with handle
- Easier to maneuver outdoors and indoors than standard Rollators with 6" casters
- Weight w/ Basket is 15 lbs
- Weight Capacity is 250 lbs

	Seat Width (in.)	seat pad depth (in.)	Seat pad width (in.)	seat height (in.)	backrest height (in.)	overall height (in.)	overall open width (in.)	weight (lb.)	weight capacity (lb.)
R-4608	18	13	13.5	21	11	33~37	24	15	250

\*\* Price & specifications are subject to change without notice.



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335