



R-4602-T

2-in-1 Luxury Aluminum Lt Weight Rollator & Transporter

- A Rollator Also is Transporter
- Brand New Design
- Foldable - Easy to Put in The Trunk
- Adjustable Handle Height w/Soft Seat
- Footplates for Easy Transport
- Loop Brake System
- 8" Flat Free Casters
- Frame Color: Metallic Burgundy or Blue
- Weight Only: 21 lbs (w/o backrest).
- Weight Capacity: 250 lbs.
- One Year Manufacturer Warranty

- Specifications
- Padded Seat: 13" x 14"
- Seat Height: 22"
- Width Between Handles: 18"
- Open Overall Width: 24"
- Overall Height: 35"~39"
- Closed Overall Width: 20"
- Closed Length: 19"

Model No.	Description
R-4602-T-BD	2-in-1 Rollator & Transporter Metallic Burgundy
R-4602-T-BL	2-in-1 Rollator & Transporter Metallic Blue

* check if current, specifications subject to change without notice



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335