



T-2700

Steel transporter, detachable arms

- Medicare code E1038
- Lightweight transporter
- Detachable desk armrests with side panels.
- Swing-away, adjustable length detachable footrests with heel loops
- Comes with seat belt included.
- 8" x 1" mag wheels with solid rubber tires
- Padded Breathable nylon upholstery
- Seat comes with strap handle on both sides for easy pull.
- Frame Color: Black

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight Capacity (lb.)
T-2700N	17	16	20	17	37	25 3/4	27	250
T-2700	19	16	20	17	37	27 3/4	27	250

* Footrests weigh 2 lbs. each.

Common Options: *Not available on all models ** Price & specifications are subject to change without notice.

• Elevating Footrests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Leg Straps



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335