



ROLLATORS

Light Weight Walking Aids

Karman lightweight series rollators are designed to provide the same comfort and reliability as our line of high-quality lightweight wheelchairs. Designed with the same ultra-lightweight high strength aluminum frame that is used in our ERGO series wheelchairs, you can expect nothing but quality from line of rollators / walkers.



Baskets or bags included on all rollators provide a convenient place to store items.



Most models come standard dual integrated braking systems for quick and controlled stops.



Most models fold as well, for easy storage and convenience.

OUR WALKERS:



R-3600 (3 wheel):

The Karman R-3600 three-wheel rollators are specifically designed to give you additional stability and comfort while walking. Because of its 3 wheel design, the R-3600 is much more maneuverable than other 4 wheel models.



R-4100 (low seat):

Weighing only 11 lbs, this rollator is a perfect match for someone looking for a more compact design.



R-4200 (pushdown brake rollator):

Lightweight rollator with a unique pushdown braking system. Perfect for those who have trouble squeezing brake levers quickly.

R-4600:

A standard and lightweight affordable rollator with dual lever braking system.



R-4602T:

What happens when you combine a lightweight rollator with a comfortable Karman wheelchair? You get the R-4602T. The best of both worlds, in one lightweight package weighing only 21 lbs.

R-4700W:

An overbuilt, extra-wide rollator with a flip-down seat. Uses a high-strength steel frame and ultra-tough brakes. Carries up to 400 lbs on the seat!



R-4800:

Much like the R-4700W, only with a lighter aluminum frame. Weight capacity of 400 lbs on the seat.

KW-100 Knee Walker:

Designed to comfortably transport a person using one knee. Features high quality design, powerful loop braking system, and a comfortable 2 piece knee pad.





Karman Healthcare
19255 San Jose Ave.
City of Industry, CA 91748

Tel: 626-581-2235
www.karmanhealthcare.com
Email: davidchao@karmanhealthcare.com