

T-2700 **Transport Class**

TRANSPORTER
DETACHABLE ARMREST

SPECIFICATIONS	T-2700N	T-2700
Seat Width	17 in.	20 in.
Seat Depth	16 in.	16 in.
Seat Height	20 in.	20 in.
Back Height	17 in.	17 in.
Overall Height	37 in.	37 in.
Overall Open Width	25¾ in.	27¾ in.
Weight w/o Riggings	27 lbs.	27 lbs.
Weight Capacity	250 lbs.	250 lbs.
HCPCS Code	E1038	E1038

(SUBJECT TO CHANGE WITHOUT NOTICE)



FIND THE **PERFECT** FIT FOR YOU...

- Detachable armrests w/ side panels
- Swing-away detachable footrests w/ heel loop
- Foldable w/ folding backrest
- 17" or 19" Seat Width of nylon
- Weight of 29 lbs. w/capacity of 250 lbs.
- Black steel frame