



Ergo Lite

Ergo Lite (S-2501) *Ultra Lightweight Ergonomic Transporter*

- Lightest Transporter On The Market
- Patented S-Style Ergonomic Seat Frame
- 6061 T-6 Aircraft-grade Aluminum Only 18 lbs. (w/ footrests)
- 3-Stage handle brake: allows light to firm grip for lock
- Built in Silver Aegis Anti-bacterial Cushion
- Folding Backrest/Folding Seat for Transporting in Vehicle or Travel
- Fixed Armrests w/Concaved Armpads
- Fixed Footrests w/Extra Wide Footplates
- Pocket Behind Backrest & Small Carry Pouch on Each Armrest
- 14" Rear Polyurethane, High Tread, Flat Free Wheels
- 6" x 1" Polyurethane Front Casters
- 16" x 17" or 18" x 17" Seat Width
- Frame Color: Pearl Silver

Features for S-Ergo Series

- Ergonomic, ultra lightweight specialized frame
- Pressure relief and anti-slippage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Concaved armpads, providing comfort & support

Only 18 lbs (w/Footrests)

***U.K. ISO 7176/19 Dynamic
Crash Test Passed***

	seat width	seat depth	armrest height	seat height	back height	overall height	overall open width	folded width	overall length	weight w/o riggings	weight capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
ErgoLite-16	16	17	8	19	17	36	23	12	39	18	220
ErgoLite-18	18	17	8	19	17	36	25	12	39	18	220

**Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam, UV Resistant Canopy, Rigidifying Push-Bar, Foldable/Adjustable Headrest*



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue, City of Industry, CA 91748-1418

Order. 1-800-805-2762 / Tel. 626-581-2235 / Fax. 626-581-2335