



KM-9020L

*aluminum high strength ultra lightweight,
super compact wheelchairs*

- Ultra light weight 6061 T-6 Aircraft-grade aluminum frame
- Medicare Code: K0004
- Folding backrests
- Fixed armrests & footrests with velcro mount leg strap
- 6" flat free front casters
- 20" flat free rear wheels with polyurethane tires
- High-grade flame retardant, breathable nylon upholstery
- Weight only: 25 lbs w/footrests
- Weight capacity: 220 lbs
- Frame color: Black, Champagne

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	Weight w/footrests (lb.)	weight capacity (lb.)
KM-9020L	18	15	19	19	37	25	25	220

FOLDED DIMENSIONS:

Overall L x W x H 30" x 13" x 26"

Model No.	Description
KM-9020L	18" seat width, aluminum, high strength ultra lightweight

Frame Color Options: **Black or Champagne**

Common Options: • Seat Belt • Pneumatic Tires • Desk Tray

* Price and specifications are subject to change without notice.

KARMA  **KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335