





Unique Features for S-Ergonomic Series

- Ergonomic, ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slippage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider arm pad, providing comfort/support

S-Ergo 105

Lightweight Ergonomic Manual Wheelchair

- High strength, weighs only 27 lbs. (w footrests)
- S-Style Ergonomic Seat
- Silver 1/4" removable Anti-Bacterial cushion
- Folding backrest/ folding seat for transporting vehicle or travel.
- Fixed armrest w/ wider concaved armpads
- Fixed Footrest
- "Tube-in Center" foot-plate, assures better side leg support
- Extending push to lock manual brakes
- Pouch for carrying small items attached to upholstery behind backrest.
- 24' Rear polyurethane, high profile, flat free wheels
- 7"x1" Polyurethane front caster
- 16x17 or 18x17 Seat width
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Black Breathable Mesh Bottom
- Chair comes fully assembled except for footrests, tools included
- Weight Capacity of 250 lbs

	seat	seat	armrest	seat	back	overall	overall	folded	overall	weight w/o	weight
	width	depth	height	height	height	height	open width	width	length	riggings	capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
S-105F16	16	17	8	19	17	35	23	12	39	27	250
S-105F18	18	17	8	19	17	35	25	12	39	27	250

*Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Foldable Canopy, Rigidifying Push-Bar, Reflecting Spoke Guard, Foldable/Adjustable Headrest, Frog Leg Suspension



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335