

**karma**

# Ergonomic

Wheelchairs

## ATX

## Ergonomic Class

# ULTRA LIGHTWEIGHT ACTIVE WHEELCHAIR

SPECIFICATIONS	ATX-14	ATX-16	ATX-18
Seat Width	14 in.	16 in.	18 in.
Seat Depth	15 in.	16 in.	18 in.
Armrest Height	23 in.	23 in.	23 in.
Seat Height	14 in.	14 in.	14 in.
Back Height	Varies	Varies	Varies
Overall Height	Varies	Varies	Varies
Overall Open Width	Varies	Varies	Varies
Folded Width	14 in.	14 in.	14 in.
Overall Length	Varies	Varies	Varies
Weight w/o Riggings	Varies	15.4	Varies
Weight Capacity	275 lbs.	275 lbs.	275 lbs.
HCPCS Code	N/A	N/A	N/A

(SUBJECT TO CHANGE WITHOUT NOTICE)



THE **ULTIMATE** STATEMENT IS ONE WITH NO COMPROMISE...

- Foot plate height unlimited increment adj.
- 15.4 LB Starting Frame Weight
- Back Cane Height and Angle Adjustable
- One piece foldable footplate
- Fixed Footrests w/ Extra Wide Footplates
- Center of Gravity Adjustable
- S-Shape Frame (patented)
- Built in Silver Aegis Anti-bacterial Cushion
- Metallic Aegis Treated High Gloss Paint
- Foot plate height unlimited increment adj.

## ErgoFlight Ergonomic Class

### ULTRA LIGHTWEIGHT AIRCRAFT GRADE T6 ALUMINUM

SPECIFICATIONS	FLIGHT-16	FLIGHT-18
Seat Width	16 in.	18 in.
Seat Depth	17 in.	17 in.
Armrest Height	8 in.	8 in.
Seat Height	18½ in.	18½ in.
Back Height	17½ in.	17½ in.
Overall Width	36 in.	36 in.
Overall Open Width	25 in.	27 in.
Folded Width	15 in.	15 in.
Overall Length	38 in.	38 in.
Weight w/o Riggings	19 lbs.	19 lbs.
Weight Capacity	220 lbs.	220 lbs.
HCPCS Code	K0004	K0004

(SUBJECT TO CHANGE WITHOUT NOTICE)



### GEOMETRY AND TECHNOLOGY TO ITS **BEST**...

- S-Style Ergonomic Seat
- 6061 T-6 Aircraft-grade Aluminum Only 19.8 lbs. (w/o footrests)
- Folding Backrest/ Folding Seat for Transporting in Vehicle or Travel
- Fixed Armrests w/ Concaved Armpads
- 20" Rear Polyurethane, High Tread, Flat Free Wheels
- 6061 T-6 Aircraft-grade Aluminum Only 14.5 lb
- Patented S-Style Ergonomic Seat Frame
- Built in Silver Aegis Anti-bacterial Cushion
- 18" x 17" Seat Width
- Frame Color: Pearl Silver



## ErgoLite Ergonomic Class

### ULTRA LIGHTWEIGHT AIRCRAFT GRADE T6 ALUMINUM

SPECIFICATIONS	LITE-16	LITE-18
Seat Width	16 in.	18 in.
Seat Depth	17 in.	17 in.
Armrest Height	8 in.	8 in.
Seat Height	19 in.	19 in.
Back Height	17 in.	17 in.
Overall Height	36 in.	36 in.
Overall Open Width	23 in.	25 in.
Folded Width	12 in.	12 in.
Overall Length	39 in.	39 in.
Weight w/o Riggings	18 lbs.	18 lbs.
Weight Capacity	220 lbs.	220 lbs.
HCPCS Code	K0004	K0004

(SUBJECT TO CHANGE WITHOUT NOTICE)



### LIGHTEST TRANSPORTER ON THE MARKET...

- S-Style Ergonomic Seat
- 6061 T-6 Aircraft-grade Aluminum Only 19.8 lbs. (w/o footrests)
- Folding Backrest/ Folding Seat for Transporting in Vehicle or Travel
- Fixed Armrests w/ Concaved Armpads
- Fixed Footrests w/ Extra Wide Footplates
- 3-Stage handle brake: allows light to firm grip for lock
- Patented S-Style Ergonomic Seat Frame
- Built in Silver Aegis Anti-bacterial Cushion
- 16" x 17" or 18" x 17" Seat Width
- Frame Color: Pearl Silver

## S-105 Ergonomic Class

### ULTRA LIGHTWEIGHT MANUAL WHEELCHAIR

SPECIFICATIONS	S-105-F16	S-105-F18
Seat Width	18 in.	18 in.
Seat Depth	17 in.	17 in.
Armrest Height	8 in.	8 in.
Seat Height	19 in.	19 in.
Back Height	17 in.	17 in.
Overall Height	36 in.	36 in.
Overall Open Width	25 in.	25 in.
Folded Width	12 in.	12 in.
Overall Length	39 in.	39 in.
Weight w/o Riggings	27 in.	27 in.
Weight Capacity	250 lbs.	250 lbs.
HCPSC Code	K0004	K0004

(SUBJECT TO CHANGE WITHOUT NOTICE)



#### ERGONOMIC, ULTRA LIGHTWEIGHT SPECIALIZED FRAME...

- Anti-Bacterial upholstery
- Folding backrest
- Fixed armrest & footrest
- "Tube-in Center" foot-plate, assures better side leg support
- 24" Rear spoke, polyurethane, high profile, flat free wheels
- Curved ( like a "S") Seating Frame
- S-shaped seating frame provides anti-slippage
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control

## S-106

## Ergonomic Class

### ULTRA LIGHTWEIGHT ADJUSTABLE BACKREST

#### SPECIFICATIONS

S-106-F18

Seat Width	18 in.
Seat Depth	17 in.
Armrest Height	8 in.
Seat Height	19 in.
Back Height	17 in.
Overall Height	36 in.
Overall Open Width	25 in.
Folded Width	12 in.
Overall Length	39 in.
Weight w/o Riggings	27 in.
Weight Capacity	250 lbs.
HCPCS Code	K0004

(SUBJECT TO CHANGE WITHOUT NOTICE)



#### 3 EASY ANGLE ADJUSTMENTS FOR THE **IDEAL** POSITION

- Adjustable angle backrest in 3 adjustments - 95°, 100°, 105°
- Fixed armrests and footrest
- High strength frame
- Lightweight (27 lbs.)
- Anti-bacterial, removable/washable upholstery
- Folding precision backrest
- "Tube-in-center" footplate assures side leg support
- 24" Rear polyurethane, high-profile flat-free wheels
- Rear anti-tippers



## S-115

## Ergonomic Class

# ULTRA LIGHTWEIGHT MANUAL WHEELCHAIR

SPECIFICATIONS	S-115-F16	S-115-F18	S-115-F20
Seat Width	16 in.	18 in.	20 in.
Seat Depth	17 in.	17 in.	17 in.
Armrest Height	8 in.	8 in.	8 in.
Seat Height	19 in.	19 in.	19 in.
Back Height	17 in.	17 in.	17 in.
Overall Height	36 in.	36 in.	36 in.
Overall Open Width	23 in.	23 in.	27 in.
Folded Width	12 in.	12 in.	12 in.
Overall Length	39 in.	39 in.	39 in.
Weight w/o Riggings	25 in.	25 in.	26 in.
Weight Capacity	250 lbs.	250 lbs.	250 lbs.
HCPCS Code	K0004	K0004	K0004

(SUBJECT TO CHANGE WITHOUT NOTICE)



## PERFECT SHAPE TO FIT THE HUMAN BODY...

- High strength, light weight
- Weighs only 25 lbs. without footrests
- Optional quick release wheels pop off with push of button
- Folding backrest
- Fixed armrests and swingaway footrests
- 24" Rear spoke flat-free wheels
- S-shaped seating frame provides anti-slippage
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control

## S-125

## Ergonomic Class

# ULTRA LIGHTWEIGHT MANUAL WHEELCHAIR

SPECIFICATIONS	S-125-F16	S-125-F18
Seat Width	16 in.	18 in.
Seat Depth	17 in.	17 in.
Armrest Height	8 in.	8 in.
Seat Height	19 in.	19 in.
Back Height	17 in.	17 in.
Overall Height	36 in.	36 in.
Overall Open Width	23 in.	23 in.
Folded Width	12 in.	12 in.
Overall Length	39 in.	39 in.
Weight w/o Riggings	25 in.	25 in.
Weight Capacity	250 lbs.	250 lbs.
HCPSC Code	K0004	K0004

(SUBJECT TO CHANGE WITHOUT NOTICE)



## PERFORMANCE, STYLE AND COMFORT...

- Fixed padded full armrests with alloy side guard
- Swing-away detachable footrests with heel loop
- Push to lock wheel breaks
- 7" Flat free polyurethane casters
- 24" Rear-Spoked, polyurethane, high tread, rear tires
- High-grade/Padded nylon flame retardant upholstery
- Seat width: 16", 18"
- Frame Color: Pearl Silver & Rose Red



## S-305

## Ergonomic Class

# ULTRA LIGHTWEIGHT MANUAL WHEELCHAIR

SPECIFICATIONS	S-305-16	S-305-18
Seat Width	16 in.	18 in.
Seat Depth	17 in.	17 in.
Armrest Height	9-11 in.	9-11 in.
Seat Height	19 in.	19 in.
Back Height	17 in.	17 in.
Overall Height	38 in.	38 in.
Overall Open Width	24 in.	26 in.
Folded Width	12 in.	12 in.
Overall Length	43 in.	43 in.
Weight w/o Riggings	29 in.	29 in.
Weight Capacity	250 lbs.	250 lbs.
HCPCS Code	K0005	K0005

(SUBJECT TO CHANGE WITHOUT NOTICE)



## ULTRA LIGHTWEIGHT **SPECIALIZED** FRAME...

- 24" Quick Release rear spoke/polyurethane/high profile/flat free wheels
- Axle Adjustable Seat Height ( 18", 19", 20")
- One piece foldable footplate
- Folding backrest, for transporting chair
- Swing in & away footrest for maximum safety while entering or exiting the chair
- S-Shape Frame (patented)
- Anti-Bacterial Upholstery/Cushion
- Frame Color: Pearl Silver & Rose Red
- Seat Width: 18"x17" OR 16" x 17"

**karma**



**Karman Healthcare**

**19255 San Jose Ave.**

**City of Industry, CA**

**Telephone Number - 626-581-2235**

**[www.karmanhealthcare.com](http://www.karmanhealthcare.com)**